



# PPEF1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA018435LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O14829
<b>Immunogen</b>	Recombinant Human Serine/threonine-protein phosphatase with EF-hands 1 protein (301-500AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Serine/threonine-protein phosphatase with EF-hands 1 (PPEF-1) (EC 3.1.3.16) (Protein phosphatase with EF calcium-binding domain) (PPEF) (Serine/threonine-protein phosphatase 7) (PP7), PPEF1, PPEF PPP7C
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	PPEF1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.