



PPEF1 Antibody, HRP conjugated

Product Code	CSB-PA018435LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14829
Immunogen	Recombinant Human Serine/threonine-protein phosphatase with EF-hands 1 protein (301-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase with EF-hands 1 (PPEF-1) (EC 3.1.3.16) (Protein phosphatase with EF calcium-binding domain) (PPEF) (Serine/threonine-protein phosphatase 7) (PP7), PPEF1, PPEF PPP7C
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPEF1