




# PPFIBP2 Antibody

<b>Product Code</b>	CSB-PA847688LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8ND30
<b>Immunogen</b>	Recombinant Human Liprin-beta-2 protein (101-400AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Liprin-beta-2 (Protein tyrosine phosphatase receptor type f polypeptide-interacting protein-binding protein 2) (PTPRF-interacting protein-binding protein 2), PPFIBP2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	PPFIBP2

## Image



Western blot  
 All lanes: PPFIBP2 antibody at 6µg/ml + HepG2 whole cell lysate  
 Secondary  
 Goat polyclonal to rabbit at 1/10000 dilution  
 Predicted band size: 99, 87, 84 kDa  
 Observed band size: 99 kDa

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.