



PPFIBP2 Antibody, HRP conjugated

Product Code	CSB-PA847688LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8ND30
Immunogen	Recombinant Human Liprin-beta-2 protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Liprin-beta-2 (Protein tyrosine phosphatase receptor type f polypeptide-interacting protein-binding protein 2) (PTPRF-interacting protein-binding protein 2), PPFIBP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPFIBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.