



PPIH Antibody, Biotin conjugated

Product Code	CSB-PA018478LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43447
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase H protein (2-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase H (PPIase H) (EC 5.2.1.8) (Rotamase H) (Small nuclear ribonucleoprotein particle-specific cyclophilin H) (CypH) (U-snRNP-associated cyclophilin SnuCyp-20) (USA-CYP), PPIH, CYP20 CYPH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPIH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.