




PPP3R1 Antibody

Product Code	CSB-PA018575LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63098
Immunogen	Recombinant Human Calcineurin subunit B type 1 protein (11-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcineurin subunit B type 1 (Protein phosphatase 2B regulatory subunit 1) (Protein phosphatase 3 regulatory subunit B alpha isoform 1), PPP3R1, CNA2 CNB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PPP3R1

Image



Western blot

All lanes: PPP3R1 antibody at 1µg/ml

Lane 1: K562 whole cell lysate

Lane 2: HeLa whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 20 kDa

Observed band size: 20 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.