



PQLC2 Antibody, Biotin conjugated

Product Code	CSB-PA744399LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZP29
Immunogen	Recombinant Human Lysosomal amino acid transporter 1 homolog protein (91-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysosomal amino acid transporter 1 homolog (PQ-loop repeat-containing protein 2), PQLC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PQLC2