





PLPBP Antibody

Product Code	CSB-PA018755LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94903
Immunogen	Recombinant Human Pyridoxal phosphate homeostasis protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyridoxal phosphate homeostasis protein (PLP homeostasis protein) (Proline synthase co-transcribed bacterial homolog protein) (Pyridoxal phosphate-binding protein), PLPBP, PROSC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLPBP

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA018755LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA018755LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.