



PLPBP Antibody, FITC conjugated

Product Code	CSB-PA018755LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94903
Immunogen	Recombinant Human Pyridoxal phosphate homeostasis protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyridoxal phosphate homeostasis protein (PLP homeostasis protein) (Proline synthase co-transcribed bacterial homolog protein) (Pyridoxal phosphate-binding protein), PLPBP, PROSC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLPBP