



Psat1 Antibody, Biotin conjugated

Product Code	CSB-PA860801HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99K85
Immunogen	Recombinant Mouse Phosphoserine aminotransferase protein (1-370AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoserine aminotransferase (PSAT) (EC 2.6.1.52) (Endometrial progesterone-induced protein) (EPIP) (Phosphohydroxythreonine aminotransferase), Psat1, Psa Psat
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Psat1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.