



PSD4 Antibody, Biotin conjugated

Product Code	CSB-PA843332LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDX1
Immunogen	Recombinant Human PH and SEC7 domain-containing protein 4 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PH and SEC7 domain-containing protein 4 (Exchange factor for ADP-ribosylation factor guanine nucleotide factor 6 B) (Exchange factor for ARF6 B) (Pleckstrin homology and SEC7 domain-containing protein 4) (Telomeric of interleukin-1 cluster protein), PSD4, EFA6B TIC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PSD4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.