




PSG11 Antibody

Product Code	CSB-PA891577LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ72
Immunogen	Recombinant Human Pregnancy-specific beta-1-glycoprotein 11 protein (7-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pregnancy-specific beta-1-glycoprotein 11 (PS-beta-G-11) (PSBG-11) (Pregnancy-specific glycoprotein 11) (Pregnancy-specific beta-1-glycoprotein 13) (PS-beta-G-13) (PSBG-13) (Pregnancy-specific glycoprotein 13), PSG11, PSG13 PSG14
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PSG11

Image



Western blot
 All lanes: PSG11 antibody at 12µg/ml + Hela whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 38, 24 kDa
 Observed band size: 38 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.