



# PSG11 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA891577LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UQ72
<b>Immunogen</b>	Recombinant Human Pregnancy-specific beta-1-glycoprotein 11 protein (7-219AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pregnancy-specific beta-1-glycoprotein 11 (PS-beta-G-11) (PSBG-11) (Pregnancy-specific glycoprotein 11) (Pregnancy-specific beta-1-glycoprotein 13) (PS-beta-G-13) (PSBG-13) (Pregnancy-specific glycoprotein 13), PSG11, PSG13 PSG14
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Immunology
<b>Target Names</b>	PSG11
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.