



PSMD7 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA018912LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P51665 |
| Immunogen | Recombinant Human 26S proteasome non-ATPase regulatory subunit 7 protein (1-324AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | 26S proteasome non-ATPase regulatory subunit 7 (26S proteasome regulatory subunit RPN8) (26S proteasome regulatory subunit S12) (Mov34 protein homolog) (Proteasome subunit p40), PSMD7, MOV34L |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | PSMD7 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |