



PSTPIP2 Antibody, Biotin conjugated

Product Code	CSB-PA884494LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H939
Immunogen	Recombinant Human Proline-serine-threonine phosphatase-interacting protein 2 protein (28-312AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proline-serine-threonine phosphatase-interacting protein 2 (PEST phosphatase-interacting protein 2), PSTPIP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PSTPIP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.