



PTPRN Antibody, Biotin conjugated

Product Code	CSB-PA613703LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16849
Immunogen	Recombinant Human Receptor-type tyrosine-protein phosphatase-like N protein (35-575AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-type tyrosine-protein phosphatase-like N (R-PTP-N) (Islet cell antigen 512) (ICA 512) (Islet cell autoantigen 3) (PTP IA-2) [Cleaved into: ICA512-N-terminal fragment (ICA512-NTF); ICA512-transmembrane fragment (ICA512-TMF); ICA512-cleaved cytosolic fragment (ICA512-CCF)], PTPRN, ICA3 ICA512
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTPRN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.