



PYCR2 Antibody, Biotin conjugated

Product Code	CSB-PA822188LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C36
Immunogen	Recombinant Human Pyrroline-5-carboxylate reductase 2 protein (2-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyrroline-5-carboxylate reductase 2 (P5C reductase 2) (P5CR 2) (EC 1.5.1.2), PYCR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.