



RAB11B Antibody, FITC conjugated

Product Code	CSB-PA614903LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15907
Immunogen	Recombinant Human Ras-related protein Rab-11B protein (2-215AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-11B (GTP-binding protein YPT3), RAB11B, YPT3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RAB11B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.