



RAB35 Antibody, Biotin conjugated

Product Code	CSB-PA614887LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15286
Immunogen	Recombinant Human Ras-related protein Rab-35 protein (1-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-35 (GTP-binding protein RAY) (Ras-related protein Rab-1C), RAB35, RAB1C RAY
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RAB35
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.