



RAB38 Antibody, Biotin conjugated

Product Code	CSB-PA019191LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57729
Immunogen	Recombinant Human Ras-related protein Rab-38 protein (1-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-38 (Melanoma antigen NY-MEL-1), RAB38
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAB38
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.