



RAB3IL1 Antibody, Biotin conjugated

Product Code	CSB-PA837424LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TBN0
Immunogen	Recombinant Human Guanine nucleotide exchange factor for Rab-3A protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide exchange factor for Rab-3A (Rab-3A-interacting-like protein 1) (Rab3A-interacting-like protein 1) (Rabin3-like 1), RAB3IL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RAB3IL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.