



# RAB3IL1 Antibody, FITC conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA837424LC01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q8TBN0  |
| <b>Immunogen</b>           | Recombinant Human Guanine nucleotide exchange factor for Rab-3A protein (1-200AA)   |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | FITC  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Guanine nucleotide exchange factor for Rab-3A (Rab-3A-interacting-like protein 1) (Rab3A-interacting-like protein 1) (Rabin3-like 1), RAB3IL1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Research Area</b>       | Neuroscience  |
| <b>Target Names</b>        | RAB3IL1   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |