



RAB40B Antibody, Biotin conjugated

Product Code	CSB-PA615534LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12829
Immunogen	Recombinant Human Ras-related protein Rab-40B protein (1-278AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-40B (SOCS box-containing protein RAR) (Protein Rar), RAB40B, SEC4L
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RAB40B