



RAB40B Antibody, Biotin conjugated

Product Code	CSB-PA615534LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12829
Immunogen	Recombinant Human Ras-related protein Rab-40B protein (1-278AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-40B (SOCS box-containing protein RAR) (Protein Rar), RAB40B, SEC4L
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RAB40B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.