



RABGGTA Antibody, Biotin conjugated

Product Code	CSB-PA852898LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92696
Immunogen	Recombinant Human Geranylgeranyl transferase type-2 subunit alpha protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Geranylgeranyl transferase type-2 subunit alpha (EC 2.5.1.60) (Geranylgeranyl transferase type II subunit alpha) (Rab geranyl-geranyltransferase subunit alpha) (Rab GG transferase alpha) (Rab GGTase alpha) (Rab geranylgeranyltransferase subunit alpha), RABGGTA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RABGGTA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.