



RALYL Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA771422LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q86SE5 |
| Immunogen | Recombinant Human RNA-binding Raly-like protein (1-291AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | RNA-binding Raly-like protein (hRALYL) (Heterogeneous nuclear ribonucleoprotein C-like 3) (hnRNP core protein C-like 3), RALYL, HNRPCL3 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | RALYL |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |