



RALYL Antibody, HRP conjugated

Product Code	CSB-PA771422LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SE5
Immunogen	Recombinant Human RNA-binding Raly-like protein (1-291AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding Raly-like protein (hRALYL) (Heterogeneous nuclear ribonucleoprotein C-like 3) (hnRNP core protein C-like 3), RALYL, HNRPCL3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RALYL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.