



RAP1A Antibody, HRP conjugated

Product Code	CSB-PA019321LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62834
Immunogen	Recombinant Human Ras-related protein Rap-1A protein (1-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rap-1A (C21KG) (G-22K) (GTP-binding protein smg p21A) (Ras-related protein Krev-1), RAP1A, KREV1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAP1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.