




RAP1B Antibody

Product Code	CSB-PA03269A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61224
Immunogen	Recombinant Human Ras-related protein Rap-1b protein (1-181AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rap-1b (GTP-binding protein smg p21B), RAP1B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAP1B

Image



Western blot

All lanes: RAP1B antibody at 0.7µg/ml

Lane 1: Mouse brain tissue

Lane 2: A549 whole cell lysate


Lane 3: A431 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 21, 16, 19, 17 kDa

Observed band size: 21 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA03269A0Rb at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.