





RAPGEF5 Antibody

Product Code	CSB-PA821676LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92565
Immunogen	Recombinant Human Rap guanine nucleotide exchange factor 5 protein (201-444AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rap guanine nucleotide exchange factor 5 (Guanine nucleotide exchange factor for Rap1) (M-Ras-regulated Rap GEF) (MR-GEF) (Related to Epac) (Repac), RAPGEF5, GFR KIAA0277 MRGEF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RAPGEF5

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA821676LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA821676LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.