



RAPGEF5 Antibody, FITC conjugated

Product Code	CSB-PA821676LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92565
Immunogen	Recombinant Human Rap guanine nucleotide exchange factor 5 protein (201-444AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rap guanine nucleotide exchange factor 5 (Guanine nucleotide exchange factor for Rap1) (M-Ras-regulated Rap GEF) (MR-GEF) (Related to Epac) (Repac), RAPGEF5, GFR KIAA0277 MRGEF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RAPGEF5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.