



RBBP9 Antibody, Biotin conjugated

Product Code	CSB-PA019394LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75884
Immunogen	Recombinant Human Putative hydrolase RBBP9 protein (1-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative hydrolase RBBP9 (EC 3.-.-.-) (B5T-overexpressed gene protein) (Protein BOG) (Retinoblastoma-binding protein 10) (RBBP-10) (Retinoblastoma-binding protein 9) (RBBP-9), RBBP9, BOG RBBP10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RBBP9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.