





RBM41 Antibody

Product Code	CSB-PA850315LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96IZ5
Immunogen	Recombinant Human RNA-binding protein 41 protein (1-413AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 41 (RNA-binding motif protein 41), RBM41
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RBM41

Image



Western blot
All lanes: RBM41 antibody at 1.5µg/ml + HeLa whole cell lysate
Secondary
Goat polyclonal to rabbit at 1/10000 dilution
Predicted band size: 48, 47, 36 kDa
Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA850315LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.