



RBM41 Antibody, Biotin conjugated

Product Code	CSB-PA850315LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96IZ5
Immunogen	Recombinant Human RNA-binding protein 41 protein (1-413AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 41 (RNA-binding motif protein 41), RBM41
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RBM41
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.