




RBP5 Antibody

Product Code	CSB-PA019484LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82980
Immunogen	Recombinant Human Retinol-binding protein 5 protein (1-135AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinol-binding protein 5 (Cellular retinol-binding protein III) (CRBP-III) (HRBPiso), RBP5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RBP5

Image



Western blot
All lanes: RBP5 antibody at 8µg/ml + Mouse liver tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 16 kDa
Observed band size: 20 kDa