



REEP2 Antibody, Biotin conjugated

Product Code	CSB-PA887109LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRK0
Immunogen	Recombinant Human Receptor expression-enhancing protein 2 protein (78-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor expression-enhancing protein 2, REEP2, C5orf19 SGC32445
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	REEP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.