



RFX6 Antibody, HRP conjugated

Product Code	CSB-PA812809LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8HWS3
Immunogen	Recombinant Human DNA-binding protein RFX6 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein RFX6 (Regulatory factor X 6) (Regulatory factor X domain-containing protein 1), RFX6, RFXDC1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RFX6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.