



Rgn Antibody, Biotin conjugated

Product Code	CSB-PA019630LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03336
Immunogen	Recombinant Rat Regucalcin protein (1-299AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regucalcin (RC) (Gluconolactonase) (GNL) (EC 3.1.1.17) (Senescence marker protein 30) (SMP-30), Rgn, Smp30
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Rgn
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.