



RGP1 Antibody

Product Code	CSB-PA838781LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92546
Immunogen	Recombinant Human RAB6A-GEF complex partner protein 2 protein (56-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAB6A-GEF complex partner protein 2 (Retrograde Golgi transport protein RGP1 homolog), RGP1, KIAA0258
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.