



RGP1 Antibody, Biotin conjugated

Product Code	CSB-PA838781LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92546
Immunogen	Recombinant Human RAB6A-GEF complex partner protein 2 protein (56-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAB6A-GEF complex partner protein 2 (Retrograde Golgi transport protein RGP1 homolog), RGP1, KIAA0258
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGP1