



RGP1 Antibody, FITC conjugated

Product Code	CSB-PA838781LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92546
Immunogen	Recombinant Human RAB6A-GEF complex partner protein 2 protein (56-350AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAB6A-GEF complex partner protein 2 (Retrograde Golgi transport protein RGP1 homolog), RGP1, KIAA0258
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.