



RHEBL1 Antibody, FITC conjugated

Product Code	CSB-PA855034LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAI7
Immunogen	Recombinant Human GTPase RhebL1 protein (1-180AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTPase RhebL1 (Ras homolog enriched in brain like-1 c) (RhebL1c) (Ras homolog enriched in brain-like protein 1) (Rheb-like protein 1) (Rheb2), RHEBL1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RHEBL1