



RIPK3 Antibody, Biotin conjugated

Product Code	CSB-PA897497LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y572
Immunogen	Recombinant Human Receptor-interacting serine/threonine-protein kinase 3 protein (1-518AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-interacting serine/threonine-protein kinase 3 (EC 2.7.11.1) (RIP-like protein kinase 3) (Receptor-interacting protein 3) (RIP-3), RIPK3, RIP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RIPK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.