



RIT2 Antibody, HRP conjugated

Product Code	CSB-PA860769LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99578
Immunogen	Recombinant Human GTP-binding protein Rit2 protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTP-binding protein Rit2 (Ras-like protein expressed in neurons) (Ras-like without CAAX protein 2), RIT2, RIN ROC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RIT2