



RMND5B Antibody, Biotin conjugated

Product Code	CSB-PA856939LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96G75
Immunogen	Recombinant Human Protein RMD5 homolog B protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein RMD5 homolog B, RMND5B
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RMND5B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.