



RNF138 Antibody, FITC conjugated

Product Code	CSB-PA851971LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVD3
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF138 protein (2-245AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF138 (EC 2.3.2.27) (Nemo-like kinase-associated RING finger protein) (NLK-associated RING finger protein) (hNARF) (RING finger protein 138) (RING-type E3 ubiquitin transferase RNF138), RNF138, NARF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RNF138
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.