



# RNF141 Antibody, Biotin conjugated

Product Code	CSB-PA019839LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVD5
Immunogen	Recombinant Human RING finger protein 141 protein (1-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RING finger protein 141 (Zinc finger protein 230), RNF141, ZNF230
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RNF141