





RNF181 Antibody

Product Code	CSB-PA019860LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0P0
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF181 protein (2-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF181 (EC 2.3.2.27) (RING finger protein 181) (RING-type E3 ubiquitin transferase RNF181), RNF181
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF181

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA019860LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA019860LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.