



RNF181 Antibody, Biotin conjugated

Product Code	CSB-PA019860LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0P0
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF181 protein (2-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF181 (EC 2.3.2.27) (RING finger protein 181) (RING-type E3 ubiquitin transferase RNF181), RNF181
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF181
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.