



RNF182 Antibody, HRP conjugated

Product Code	CSB-PA822714LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6D2
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF182 protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF182 (EC 2.3.2.27) (RING finger protein 182) (RING-type E3 ubiquitin transferase RNF182), RNF182
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF182
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.